


# TL-2800 ASX

## @ the Edge of Excellence

Timberline Instruments - ASX Series


Fast - Accurate - Precise

### TL-2800 ASX NH<sub>3</sub>/Nitrate Analyzer



### EPA Approved Method for Clean Water Act Applications

"Determination of Inorganic Ammonia by Continuous Flow Gas Diffusion and Conductivity Cell Analysis" Timberline Method Ammonia-001





Gas Diffusion and Conductivity Measurement Method

The Model TL-2800 ASX Ammonia and Nitrate & Nitrite Analyzer provides automated analysis of ammonia and nitrate/nitrite in drinking water, surface water, wastewater and other aqueous samples including agricultural extracts from soil, plants and wastes.

Our patented gas diffusion & conductivity principle provides results **unaffected by turbidity or color**. Simplicity and ease of use makes the TL-2800 extremely popular.

Standard analysis ranges from 50 ppb to 500 ppm. **Ultra-sensitive detection** (<10ppb) is easily achieved with a high sensitivity membrane.

The TL-2800A ASX produces total NH<sub>3</sub>/NH<sub>4</sub><sup>+</sup> results and when paired with the TL-2800N ASX which incorporates a reduction catalyst, the analysis provides total inorganic nitrogen. Nitrate/nitrite are determined by subtraction.



1880 S. Flatiron Ct, Unit G  
Boulder, Co 80301  
800-777-5996  
info@timberlineinstruments.com

## Why choose Timberline TL-2800 ASX?

Integral Membrane Separation & Conductivity Detection

Results are unaffected by turbidity or color of the sample. Significantly reduces interferences.

Low Maintenance

Samples are introduced via peristaltic pump with timed injection. The absence of sample injection valves reduces the need for filtration and prevents costly repairs and maintenance issues associated with premature valve failure.

High Sample Throughput

A 1.5-2.5 minute analysis time allows for rapid determination, reducing the total cost per sample.

Minimal Reagent Requirements  
Minimal Disposal Issues

Only three reagents are required for operation.

- Dilute Borate Buffer
- Dilute Sodium Hydroxide
- Deionized Water

Activated zinc cartridges are used for nitrate analysis.

Familiar Data Acquisition  
System Hardware

Data acquisition system software runs on Microsoft® Windows operating system, allowing for easy connection to networks, printers or external USB drives.

Intuitive Software

The Microsoft® Windows based data acquisition software provides complete automation of the instrument and sampling system, while the intuitive graphical user interface makes it very easy to learn and operate.

High Sensitivity

Ultra sensitive detection down to less than 10 ppb.

## Specifications:

Dimensions/Instrument (L x W x H)

18" x 9.5" x 7.5" (460mm x 240mm x 190mm)

Weight

9 lbs (4 Kg)

Electrical Power

120/240VAC @ 50/60Hz - 150W

Analysis Time

One sample every 1.5-2.5 minutes

Dynamic Range

Standard Sensitivity - 50 ppb to 500 ppm  
High Sensitivity - 10 ppb to 10 ppm

Detection Limit

Standard Sensitivity - 14 ppb  
High Sensitivity - 2 ppb

Precision and Accuracy

(Ammonia) ± 5% of the amount present  
(Nitrate/Nitrite) ± 10% of the amount present

Calibration Frequency

Once per day

