

AUTOSAMPLERS

@ the Edge of Excellence

Timberline Instruments - ASX Series

Automating Your Analyses

Up to 240 Sample Tubes



Compatible Autosamplers from CETAC

CETAC ASX-280

- 120 Position Autosampler
- Holds 120 x 14ml sample tubes and 10 x 50ml tubes for standards



CETAC ASX-560

- 240 Position Autosampler
- Holds 240 x 14ml sample tubes and 10 x 50ml tubes for standards



 **TIMBERLINE**
INSTRUMENTS

1880 S. Flatiron Ct, Unit G
Boulder, Co 80301
800-777-5996

info@timberlineinstruments.com